

### PXS5-T Series Zoom Microscope

PXS5 series zoom microscope adopts streamline integrated mechanical design, Easy operation, long working distance, clear resolved image, enclosure integrated. The microscope is used extensively world medical research and health care, biology and botany research, and agriculture, as well as in electronic component manufacturing . Specially suited for LED, PCB inspection, pressing plating and electronic component inspection.

Configuration: Pillar base with reflected and transmitted LED illumination

#### **Main Specification:.**

Binocular&Trinocular Head Angle of viewing :45°

Zoom ratio :1:6.4

Eyepiece: WF10X/20mm

Objective :0.7x-4.5x

Magnification:7x-45x

Interpupillary distance adjustment range :55-75mm

Working distance:100mm

